

# AIRMODUS

Do you need to detect the particle formation?

Do you need to monitor the smallest particles in your process - online?

Or maybe resolve if your environment contains precursors for nanoparticles?



installation

✓ **lab / stationary**

compliance

✓ **ACTRIS / other**

setup

✓ **PSM + CPC**

size distribution

✓ **1 - 4 nm**

how to purchase

✓ **buy / rent**



Airmodus A11 includes unique Particle Size Magnifier which grows the particles for a CPC to count them. For A11 you do not lose sensitivity and you measure both, neutral and charged clusters! With A11 you can detect particle concentrations from 1 to 300 000 #/cc!

Let's plan a try-out based on your measurement needs!



**Enviro**

Technology Services Ltd

Tel: +44 (0) 1453 733200

[sales@et.co.uk](mailto:sales@et.co.uk)

[www.et.co.uk](http://www.et.co.uk)

Kingfisher Business Park, London Road, Stroud, Gloucestershire, GL5 2BY, UK